
# ASTOR







## PERFECT INTEGRATION

The characteristic features of LAP laser systems are sophisticated technology, quality and design for more than 30 years. This level of excellence has made us the global market leader for patient alignment in radiotherapy.

In radiotherapy precise, reproducible and fast positioning of patients are required for successful treatments and high patient throughput.

After the initial patient marking process at CT or MRI the patient will be repositioned on the LINAC using cross hair and sagittal line lasers. The fixed lasers are aligned to the isocenter of the LINAC and enable isocentric positioning of the patient.

Integrated in a superior aluminium housing our ASTOR room laser combines advanced optics, high precision mechanics and state of the art electronics to ensure the highest level of stability and reliability.

Offering various configurations the ASTOR lasers are perfectly suited to meet all room requirements.

#### LAP – WE KNOW PATIENT ALIGNMENT

- Since 1984
- Global market leader
- In-house hardware and software development
- Scientific collaborations
- Worldwide service network
- Certified in accordance to ISO 9001 and ISO 13485
- Made in Germany





#### COMPACT DESIGN

ULTRA-FINE PRECISE LINES

**DISTORTION-FREE WINDOWS** 

## ASTOR - UNIQUE TECHNOLOGY

#### CHOOSE YOUR LASER COLOR

All ASTOR room lasers come equipped with long-life red or green laser diodes. Depending on your personal preferences, ambient room light or room design you can choose your laser color.




#### PRECISE

LAP ULTRALINE is the result of advanced mechanical components and unique optoelectronics used to generate and align laser lines for medical applications. The ultra-fine and long lines meet the highest quality requirements for linearity and brightness distribution.



#### DISTORTION-FREE

LAP ASTOR lasers are fitted with unbreakable, specially flattened glass windows. This minimizes scattering and guarantees ultra-fine lines at all transition angles.





#### EASY TILT ADJUSTMENT


The ASTOR laser's horizontal and vertical tilt can be easily adjust without removing the cover. Translation, rotation and focus are adjusted after removing the cover via the convenient quick lock mechanism.

- 1 Tilt vertical plane
- 2 Tilt horizontal plane



#### MOUNTING VERSIONS

For patient alignment each LINAC room is equipped with at least three room lasers. One laser projects a sagittal line and two cross hair lasers project the transverse and horizontal lines. You are free to select from four mounting versions to perfectly match your existing room situation.








#### **DID YOU KNOW?**

To avoid shadowing effects from the LINAC gantry a fourth ASTOR cross hair can be installed overhead.

# ASTOR - HIGH PRECISION TOOLS FOR PATIENT POSITIONING







plates and the false ceiling or wall recess.

#### TILT BRACKET

LAP ASTOR lasers are equipped with adjustable tilting brackets allowing for rotation up to 45°. The ASTOR's compact design makes it an ideal solution for tight mounting locations.

# FREE STANDING COLUMN

For locations where wall mounts are not possible the LAP ASTOR laser can be equipped with an optional free standing column.



| TECHNICAL DATA                       |                              |
|--------------------------------------|------------------------------|
| Laser color                          | red (638 nm), green (520 nm) |
| Laser class                          | 2                            |
| Line width up to 4 m distance        | < 1 mm                       |
| Line length at 3 m distance          | > 3 m                        |
| Power supply                         | 100 240 V AC                 |
| Internal voltage                     | 24 V DC                      |
| Operating Temperature                | 15 30 °C                     |
| Dimensions ( $H \times W \times D$ ) | 182 × 86 × 90 mm             |
| Weight                               | 1.4 kg                       |









LAP ASTOR® is a registered trademark of LAP GmbH Laser Applikationen. Further designations of products or services may be registered trademarks of LAP  $\mathsf{GmbH}$ or other organizations; their use by third parties may infringe the rights of the respective owners.

#### LAP GmbH Laser Applikationen

Zeppelinstrasse 23 21337 Lueneburg Germany

Phone +49 4131 9511-95 +49 4131 9511-96

Email info@lap-laser.com

#### LAP of America, LLC

161 Commerce Rd., Suite 3 Boynton Beach, FL 33426 USA

Phone +1 561 416-9250 +1 561 416-9263

Email america@lap-laser.com

#### LAP GmbH Laser Applikationen

Представительство в Москве

1, Казачий переулок 7 119017 Москва Российская Федерация

+7 495 7304043 Факс +7 495 7304044

Email info-russia.hc@lap-laser.com

#### **LAP Laser Applications** Asia Pacific Pte. Ltd.

750A Chai Chee Road #07-07 Viva Business Park Singapore 469001

Phone +65 6536 9990 +65 6533 6697

Email info-asia.hc@lap-laser.com

#### **LAP Laser Applications** China Co. Ltd.

East Unit , 4F Building # 10 LuJiaZui Software Park No. 61 Lane 91 EShan Road Shanghai 200127

China

Phone +86 21 5047-8881 Fax +86 21 5047-8887 Email info-cn.hc@lap-laser.com

