

Accessories

APOLLO MR3T

Included



Remote control

- All degrees of freedom can be adjusted via remote control: shift, tilt, rotation and focus



Central power supply

- Enables central switching of up to 4 devices in one room
- Including one power supply, T-Connectors as needed, cables as needed
- Cable length from central switch to device: 15 m
- Including cable package for connection to RF-filter

Optional



Free standing column

- For locations where wall mounts are not possible
- Length options: 900–1.400 mm or 1.100–1.600 mm



AQUARIUS Phantom

- MRI phantom for external laser alignment
- 3 leveling aids
- MR-safe

All accessories are MR conditional up to 3 Tesla

LAP GmbH Laser Applikationen
 Zeppelinstr. 23
 21337 Lüneburg
 Germany

P +49 4131 95 11-95
 E info@lap-laser.com

made
 in
 Germany