

Accessories

DORADO

Included



WILKE Phantom ■ For regular laser QA

Optional



Mounting set for suspended ceiling ■ Stable bracket and mounting plate with multiple screw options for large distances between ceiling and device
■ Length up to 1.000 mm



DORADO protect ■ Protective mount to protect DORADO from external impacts, e.g. caused by patient beds, wheelchairs or positioning aids
■ Length options: 500 mm or 1.820 mm



Base for double floor ■ Mount for rooms with framed or double floors
■ Length up to 500 mm



Swivel mount for DORADO ■ Bracket for tilting the ceiling device
■ Swivel angle up to 55°



Free standing column* ■ For locations where wallmounts for ASTOR are not possible
■ Length options (isocentre height): 900–1.400 mm or 1.100–1.600 mm



2nd swivel mount* ■ Bracket to enable 2nd tilting axis for ASTOR
■ i.e. for mounting sagittal lasers at diagonally placed machines
■ Swivel angle up to 45°

* Only available for DORADO 1

LAP GmbH Laser Applikationen
Zeppelinstr. 23
21337 Lüneburg
Germany

P +49 4131 95 11-95
E info@lap-laser.com

made
in
Germany