

### **Data Sheet**

## **APOLLO**

# Remotely adjustable fixed laser system for patient positioning at LINAC

#### **Features**

- Precise laser line due to a max. width of 0.5 mm, a line straightness of ±0.1 mm and distortion-free glass with anti-reflective coating.
- Stable thanks to a high-quality aluminum housing which reliably protects the fine mechanics and optical components inside from vibration, damage and radiation exposure.
- Flexible for your room situation: Install your lasers safely, practically and aesthetically with the optional accessories.
- Comfortable laser adjustment via remote control:
  All lasers can easily be selected and adjusted in focus, rotation, tilt or shift.

#### **Configurations**

- Color (each laser): red, green or blue
- Line shape (each laser): line, tilted line 90° or crosshair
- Preconfigured sets with 3 or 4 lasers available

#### **Mounting options**









wall-ceiling-wall

post-ceiling-post

post-ceiling-wall wall-

wall-ceiling-post

#### Laser

Dimensions (L × W × H)	237 × 110 × 105 mm 9.3 × 4.3 × 4.1 "
Weight	2.6 kg
Power supply	100-240 V AC, 50-60 Hz
Power supply (internal)	24 V DC
Laser color	red (638 nm), green (520 nm), blue (450 nm)
Focusable range	1–4 m
Line length at 3 m distance	>3 m
Line width up to 4 m distance	<0.5 mm (blue), <1 mm (red, green)
Swivel angle	45°
Ambient conditions	35-80 % rel. humidity, non-condensing
Operating temperature	15-30 °C
International protection rating	IP 20
Approved	A CE conformity assessment according to the European Medical Device Regulation MDR – (EU) 2017/745 was carried out successfully. The device is listed in the FDA Medical Devices Database and legally marketed in the USA.
Laser class	2
Laser power	<1 mW

#### **Remote control**

Dimensions (L × W × H)	181 × 49 × 21 mm 7.1 × 1.9 × 0.8 ″
Weight	130 g
Power supply	3 V DC (2 batteries type AAA/ LR 03/Micro, 1.5 V)
International protection rating	IP 40

#### **Scope of delivery**

- Laser system as configured
- Adjustable tilting brackets
- Plug-in power supply
- Remote control for laser adjustment
- 12 months standard warranty

LAP GmbH Laser Applikationen Zeppelinstr. 23 21337 Lüneburg Germany

P +49 4131 95 11-95

E info@lap-laser.com

#### **Accessories**

Optional accessories available, please refer to acessory overview:







