

# Accessories

## APOLLO

### Included



#### Remote control

- All degrees of freedom can be adjusted via remote control: shift, tilt, rotation and focus

### Optional



#### Mounting set for suspended ceiling

- Stable bracket for large distances between mounting plates and false ceiling or wall recess
- Length up to 1.500 mm



#### Free standing column

- For locations where wall mounts are not possible
- Length options (isocentre height): 900–1.400 mm or 1.100–1.600 mm



#### 2<sup>nd</sup> swivel mount

- Bracket to enable 2<sup>nd</sup> tilting axis
- i.e. for mounting sagittal lasers at diagonally placed machines
- Swivel angle up to 45°



#### Base for double floor

- Mount for rooms with framed or double floors
- Length up to 500 mm



#### Central power supply

- Enables central switching of up to 4 devices in one room
- Including one power supply, T-Connectors as needed, cables as needed
- Cable length from central switch to device: 15 m



#### WILKE Phantom

- For regular laser QA

LAP GmbH Laser Applikationen  
 Zeppelinstr. 23  
 21337 Lüneburg  
 Germany

P +49 4131 95 11-95  
 E info@lap-laser.com

