

## Accessories

# APOLLO MR3T

### Included



#### Remote control

- Comfortable laser adjustment via remote control:  
All lasers can be easily selected and adjusted in focus, rotation, tilt or shift.
- With hand strap for MR conditional handling



#### Central power supply

- Enables central switching of up to 4 devices in one room
- Including one power supply, T-Connectors as needed, cables as needed
- Cable length from central switch to device: 15 m
- Including cable package for connection to RF-filter

### Optional



#### Free standing column

- For locations where wall mounts are not possible
- Length options (isocentre height): 900–1.400 mm or 1.100–1.600 mm



#### AQUARIUS Phantom

- MRI phantom for external laser alignment
- 3 leveling aids
- MR-safe

All accessories are MR conditional up to 3 Tesla

LAP GmbH Laser Applikationen  
 Zeppelinstr. 23  
 21337 Lüneburg  
 Germany

P +49 4131 95 11-95  
 E info@lap-laser.com

made  
 in  
 Germany