

Accessories

DORADOnova MR3T

Included



Remote control

- Comfortable laser adjustment via remote control:
All lasers can be easily selected and adjusted in focus, rotation, tilt or shift.
- With hand strap for MR conditional handling

Optional



AQUARIUS Phantom

- MRI phantom for external laser alignment
- 3 leveling aids



All accessories are MR conditional up to 3 Tesla

LAP GmbH Laser Applikationen
Zeppelinstr. 23
21337 Lüneburg
Germany

P +49 4131 95 11-95
E info@lap-laser.com

